



Amplification and Illumina Sequencing of the V4 region of the 18S rRNA gene

1.0 Introduction

The protocol detailed here is designed to amplify the V4 region of the 18S rRNA gene for paired-end 18S community sequencing on the Illumina MiSeq platform. This protocol is based on Illumina's 16S Metagenomic Sequencing Library Preparation guide and the protocol used by Ocean Sampling Day, modified to amplify the target and add indexed adapter sequences in a single PCR step.

2.0 Amplification of the V4 region of the 18S rRNA gene

2.1 Primers for amplification of V4 region of the 18S rRNA gene

Amplification primers

See Appendix 1 for full list of primer sequences.

Forward primer

Field number (space-delimited), description:

5' Illumina adapter

Nextera XT i5 index sequence

Illumina forward overhang sequence

18S V4 forward

AATGATACGGCGACCACCGAGATCTACAC XXXXXXXX TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG CCAGCASCYCGGTAATCC

Reverse primer

Field number (space-delimited), description:

Reverse complement of 3' Illumina adapter

Nextera XT i7 index sequence

Illumina reverse overhang sequence

18S V4 reverse

CAAGCAGAAGACGGCATACGAGAT XXXXXXXX GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAG ACTTTCGTTCTTGATYRATGA

2.2 Master mix for amplification of V4 region of the 18S rRNA gene

Component	Volume 1 rxn
KAPA HiFi Hot Start Readymix (2x) ^(a)	12.5
H ₂ O	9
Forward primer (10 µm)	1.25
Reverse primer (10 µm)	1.25
Template	1
Total Volume	25

Kit code KK2601 or KK2602

2.3 Thermocycler Conditions for amplification of V4 region of the 18S rRNA gene (96 well thermocyclers)

	Temperature	Time (mm:ss)
Activation	98°C	0:30
Amplification (10 cycles)	98°C	0:10
	44°C	0:30
	72°C	0:15
Amplification (20 cycles)	98°C	0:10
	62°C	0:30
	72°C	0:15
Final Extension	72°C	7:00
HOLD	4°C	∞

2.4 Process

- 2.4.1 Use undiluted DNA as a first attempt, and 1:10 diluted for repeats/failed reactions
- 2.4.2 Amplify samples with conditions outlined above.
- 2.4.3 Run amplicons on an agarose gel. Expected band size for 18S-V4 is approximately 536bp.
- 2.4.4 Clean and normalize the PCR products using SequelPrep Normalization plates according to manufacturer's instructions (Invitrogen cat no. A10510-01)
- 2.4.5 Pool equal volumes of each normalized amplicon.
- 2.4.6 Perform QC on pool using Qubit (concentration) and Tapestation (size) and calculate molarity of pool.

3.0 Sequencing of V4 region of the 18S rRNA gene

3.1 Sequencing Setup

- 3.2 Dilute pool prepared in **step 2.4.6** to **4nM**.
- 3.3 Denature according to Illumina protocol. See *Preparing Libraries for Sequencing on the MiSeq (part #15039740)*.
- 3.4 Prepare MiSeq Reagent Cartridge (v2 500-cycles). See *MiSeq Reagent Kit v2 - Reagent Preparation Guide (part # 15034097)*.
- 3.5 Load 600 µl of library pool into the MiSeq reagent cartridge in designated reservoir
- 3.6 Prepare sample sheet to include the appropriate index sequences.
- 3.7 Start sequencing run following *MiSeq System User Guide (part # 15027617)*.

References

16S Metagenomic Sequencing Library Preparation (Illumina Part # 15044223 Rev. B) available here:

http://www.illumina.com/content/dam/illumina-support/documents/documentation/chemistry_documentation/16s/16s-metagenomic-library-prep-guide-15044223-b.pdf

LifeWatch Italy Ocean Sampling Day 2014 Protocol – available here: http://mb3is.megx.net/osd-files/download?path=/2014/protocols&files=OSD2014_protocol_B_18S_V4andV9_Sequencing_LifeWath_MoBiLab_BA_RI.pdf

Appendix 1. Primer Sequences

Primer Name	Primer Sequence
18S-V4f_S502	AATGATACGGCGACCACCGAGATCTACACCTCTATTCTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S503	AATGATACGGCGACCACCGAGATCTACACTATCCTCTTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S505	AATGATACGGCGACCACCGAGATCTACACGTAAGGAGTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S506	AATGATACGGCGACCACCGAGATCTACACACTGCATATCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S507	AATGATACGGCGACCACCGAGATCTACACAAGGAGTATCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S508	AATGATACGGCGACCACCGAGATCTACACCTAAGCCTTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S510	AATGATACGGCGACCACCGAGATCTACACCGTCTAATTCTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S511	AATGATACGGCGACCACCGAGATCTACACTCTCTCCGTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S513	AATGATACGGCGACCACCGAGATCTACACTCGACTAGTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S515	AATGATACGGCGACCACCGAGATCTACACTTCTAGCTTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S516	AATGATACGGCGACCACCGAGATCTACACCTAGAGTTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S517	AATGATACGGCGACCACCGAGATCTACACGCGTAAGATCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S518	AATGATACGGCGACCACCGAGATCTACACCTATTAAGTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S520	AATGATACGGCGACCACCGAGATCTACACAAGGCTATTCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S521	AATGATACGGCGACCACCGAGATCTACACGAGCCTTATCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4f_S522	AATGATACGGCGACCACCGAGATCTACACTTATGCGATCGTCGGCAGCGTCAGATGTGTATAAGAGACAGCCAGCASCYCGGTAATTCC
18S-V4Lr_N701	CAAGCAGAAGACGGCATAACGAGATTCGCCTTAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N702	CAAGCAGAAGACGGCATAACGAGATCTAGTACGGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N703	CAAGCAGAAGACGGCATAACGAGATTTCTGCCTGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N704	CAAGCAGAAGACGGCATAACGAGATGCTCAGGAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N705	CAAGCAGAAGACGGCATAACGAGATAGGAGTCCGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N706	CAAGCAGAAGACGGCATAACGAGATCATGCCTAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N707	CAAGCAGAAGACGGCATAACGAGATGTAGAGAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N710	CAAGCAGAAGACGGCATAACGAGATCAGCCTCGGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N711	CAAGCAGAAGACGGCATAACGAGATTGCCTCTGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N712	CAAGCAGAAGACGGCATAACGAGATTCTCTACGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N714	CAAGCAGAAGACGGCATAACGAGATTACTACGCTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N715	CAAGCAGAAGACGGCATAACGAGATCCTGAGATGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N716	CAAGCAGAAGACGGCATAACGAGATTAGCGAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N718	CAAGCAGAAGACGGCATAACGAGATGTAGCTCCGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N719	CAAGCAGAAGACGGCATAACGAGATTACTACGCTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N720	CAAGCAGAAGACGGCATAACGAGATAGGCTCCGGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N721	CAAGCAGAAGACGGCATAACGAGATGCAGCGTAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N722	CAAGCAGAAGACGGCATAACGAGATCTGCGCATGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N723	CAAGCAGAAGACGGCATAACGAGATGAGCGTAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N724	CAAGCAGAAGACGGCATAACGAGATCGCTCAGTGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N726	CAAGCAGAAGACGGCATAACGAGATGTCTTAGGGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N727	CAAGCAGAAGACGGCATAACGAGATACTGATCGGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N728	CAAGCAGAAGACGGCATAACGAGATTAGCTGCAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA
18S-V4Lr_N729	CAAGCAGAAGACGGCATAACGAGATGACGTCGAGTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGACTTTCTGTTCTTGATYRATGA