

## DNA sequencing of metranscriptomics samples on the HiSeq 2500

### 1. Introduction

This protocol describes the procedure for preparing libraries from metatranscriptomic samples to be sequenced on the Illumina HiSeq 2500 platform. The library preparation follows the Ribo-Zero rRNA Removal Kit Reference Guide (Document # 15066012 v02, August 2016) and TruSeq Stranded mRNA Sample Preparation Guide (Document # 15031047 Rev. E, October 2013) with some protocol variations.

### 2. Sample QC

- 2.1. Verify the quality of the extracted RNA using the Bioanalyzer RNA 6000 Nano or RNA 6000 Pico Kit.

### 3. Library Preparation

*Note:* The Ribo-Zero rRNA Removal kits support rRNA depletion from 1-5 $\mu$ g of total RNA. Probes from the Bacteria and Plant Ribo-Zero kits can be combined (50:50 ratio) to simultaneously deplete both bacterial and plant rRNA. In this instance the inputs of RNA should be halved (i.e. 500 ng -2.5 ng of total RNA). The inputs for the sediment RNA sample may be much lower than these recommended inputs.

- 3.1. Bring  $\leq 1$   $\mu$ g of total RNA to a final volume of 28 $\mu$ l with RNase-free water.
- 3.2. Perform rRNA depletion following the Ribo-Zero Kit Reference guide
- 3.3. Perform clean-up of rRNA depleted supernatant using Agencourt RNAClean XP Kit according to the Ribo-Zero Kit Reference Guide. Elute samples in 8.5  $\mu$ l of RNase-free water
- 3.4. Perform TruSeq Stranded mRNA-seq library preparation according to the TruSeq Stranded mRNA Sample Preparation Guide, skipping the purification of poly(A) RNAs.  
Add 13  $\mu$ l of *Fragment, Prime and Finish Mix* to 5  $\mu$ l of rRNA depleted rRNA sample and follow the library prep protocol from the *Incubate RFP* step (page 20)
- 3.5. At the Enrich DNA Fragments step (page 38-41) perform 13-15 cycles of PCR  
*For inputs of 1  $\mu$ g use 13 cycles, for inputs < 200 ng use 15 cycles*

### 4. Library QC

- 4.1. Assess the size of the RNA-seq libraries via electrophoresis using the Agilent TapeStation TapeScreen DNA 1000 Assay or similar (Agilent Bioanalyzer, Perkin Elmer LabChip GXII)
- 4.2. Quantify libraries using qPCR (KAPA Library Quantification Kits for Illumina or similar)
- 4.3. Normalize the libraries and pool libraries for sequence following the library prep protocol (*Normalize and Pool Libraries*, page 44).

### 5. Sequencing

- 5.1. Denature the 2nM pool according to the Illumina's HiSeq and GAIIx Systems - Denature and Dilute Libraries Guide (Document # 15050107 v03, November 2016).
- 5.2. Following Illumina's cBot System Guide (Document # 15006165 v02 January 2016), perform clustering of HiSeq PE Cluster Kit v4 (PE-401-4001)
- 5.3. Following Illumina's HiSeq 2500 System Guide System Guide (Document # 15035786 v01 October 2015), prepare the HiSeq SBS v4 Reagents for a 2x100 bp PE or 2x125 bp PE sequencing run.
- 5.4. Start the HiSeq 2500 Rapid sequencing run according to the system guide referenced above.
- 5.5. After sequencing is complete, perform the basecalling and demultiplexing using bcl2fastq Conversion Software.